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
The screenshot shows the BioEdit Sequence Alignment Editor interface. The main window displays three DNA sequences (Pes 1, Pes 2, Pes 3) aligned. The sequences are:

```
Pes 1 dle F GCTCCCCGCGCGCCCCCGTGGTCCCGCCCCCGGGCCCCCTCCTGGCCACTCACGTGCCGACGCCGCCAGGGGCGCGGGCCGGATGCGCGGGCGGGCATTTGGCTGCGGGCTGGGGTG
Pes 2 dle F CGCGGGCCCCCGTGGTCCCGCCCCCGGGCCCCCTCCTGGCCACTCACGTGCCGACGCCGCCAGGGGCGCGGGCCGGATGCGCGGGCGGGCATTTGGCTGCGGGCTGGGGTG
Pes 3 dle F CGCGGGCCCCCGTGGTCCCGCCCCCGGGCCCCCTCCTGGCCACTCACGTGCCGACGCCGCCAGGGGCGCGGGCCGGATGCGCGGGCGGGCATTTGGCTGCGGGCTGGGGTG
```

A dialog box titled "BioEdit Sequence Alignment Editor" is open in the foreground, displaying the following text:

A log file is being generated.
It will be located in C:\BioEdit\Temp.
The log file for this session will be called BioEdit_clustalw_012021_155433.log.
The temporary input file for this run will be called ~fasta_012021_155433.tmp.
The temporary output file will be called ~align_012021_155433.tmp.
Run parameters will be in a file called clusParm_012021_155433.txt.
All of these files will be located in C:\BioEdit\Temp
Note: These files will not be deleted automatically.
You may want to clean them out of your swap directory periodically.
NOTE THE FOLLOWING IF THIS IS A TRANSLATED NUCLEOTIDE ALIGNMENT:
Any stops (***) will be replaced by 'X' for clustal.
They will revert to stops (***) upon re-toggling the translation.

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




Reverse complement SEQ R: 5' GCGCGGTGACTCAGCACT 3'

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SEQ F: 5' GTTCTGGGC **GTTTAGAGC** 3'


Reverse complement SEQ R: 5' **GCGCG**GTGACTCAGCACT 3'

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The screenshot shows the NCBI search results for the query 'Canis familiaris'. The search was performed on the NCBI website, and the results are displayed in a structured format. The main result is the 'ROS_Cfam_1.0' genome assembly, which is a reference sequence for the dog genome. The assembly was generated by The Roslin Institute in September 2020 and is available in the RefSeq database (GCF_014441545.1). The search results also include a section for 'Assembly statistics' and two columns of related data: 'Literature' and 'Genes'. The 'Literature' column lists various databases and the number of records found in each, while the 'Genes' column lists various databases and the number of records found in each.

Search NCBI

Results found in 27 databases (1 error)

GENOME ASSEMBLY Was this helpful?

ROS_Cfam_1.0
Canis lupus familiaris (dog)
The Roslin Institute (September 2020)
RefSeq GCF_014441545.1

[Genome Browser](#) [BLAST](#) [Get data](#)

Assembly statistics +

Literature	
Bookshelf	4,210
MeSH	23
NLM Catalog	1,564
PubMed	347,753
PubMed Central	186,112

Genes	
Gene	50,917
GEO DataSets	14,784
GEO Profiles	205,497
HomoloGene	17,434
PopSet	536

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NIH U.S. National Library of Medicine NCBI National Center for Biotechnology Information

BLAST » blastn suite

Canis lupus familiaris (dog) Nucleotide BLAST

blastn blasto blastx tblastn tblastx

BLASTN programs search nucleotide databases using a nucleotide query. [more...](#)

Enter Query Sequence

Enter accession number(s), gi(s), or FASTA sequence(s) [Clear](#) Query subrange [From](#) [To](#)

```
>Consensus dle F i R
GTTTATAGACCGCTGGGCGAGGAGCCGCGCCGCGGGGCGCGAGGCTCCCGCGCGGCC
CCGTGGTCCCGCCCGCCGCTCTGCGCACTACAGTGGCGAGCGCCCGAGGGCGCGGG
CCGGATGCGGGGCGCGATGGCTGCGGGCTGGGGTGGGGGCGCGCGCTCCGGGCTA
TAAGAGCGCGGGGCGCGCCGCTGGGTTGGCTGGCGTCTGCTGCGGGCTCTCCGCTGC
CCTCGGACTGCTGAGCCGTCGGAGCCCGCCGAGGAGCGCGCCCGCCCTCGCGAGTCAIGGA
```

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Job Title
Enter a descriptive title for your BLAST search

Choose Search Set

Database

Exclude Models (XM/XP)

Optional Entrez Query
Optional Enter an Entrez query to limit search

Program Selection

Optimize for

- Highly similar sequences (megablast)
- More dissimilar sequences (discontiguous megablast)
- Somewhat similar sequences (blastn)

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U.S. National Library of Medicine
National Center for Biotechnology Information

BLAST » blastn suite » results for RID-0K308KH7013

Home Recent Results Saved Strategies Help

[← Edit Search](#) Save Search Search Summary ▾

How to read this report? BLAST Help Videos Back to Traditional Results Page

Job Title Consensus dle F I R

RID 0K308KH7013 Search expires on 01-22 21:13 pm [Download All](#) ▾

Program BLASTN [?](#) [Citation](#) ▾

Database Genome (ROS_Cfam_1.0 reference, Annotation Release 106)
[See details](#) ▾

Query ID lcl|Query_27197

Description Consensus dle F I R

Molecule type dna

Query Length 570

Other reports [Distance tree of results](#) [MSA viewer](#) [?](#)

Descriptions Graphic Summary Alignments Taxonomy

Sequences producing significant alignments Download ▾ [New](#) Select columns ▾ Show 100 ▾ [?](#)

select all 1 sequences selected [GenBank](#) [Graphics](#) [Distance tree of results](#) [New](#) [MSA Viewer](#)

	Description	Scientific Name	Max Score	Total Score	Query Cover	E value	Per. Ident	Acc. Len	Accession
<input checked="" type="checkbox"/>	Canis lupus familiaris isolate SID07034 breed Labrador retriever chromosome 31 ROS_Cfam_1.0	Canis lupus familiaris	1048	1048	100%	0.0	99.82%	39518933	NC_051835.1

Filter Results

Organism only top 20 will appear exclude

Type common name, binomial, taxid or group name

[+ Add organism](#)

Percent Identity to **E value** to **Query Coverage** to

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Primer3 (v. 0.4.0) Pick primers from a DNA sequence.

There is a newer version of Primer3 available at <http://primer3.ut.ee>

Paste source sequence below (5'>3', string of ACGTNaactn -- other letters treated as N -- numbers and blanks ignored). FASTA format ok. Please N-out undesirable sequence (vector, ALUs, LINES, etc.) or use a [Mispriming Library \(repeat library\)](#): NONE

```
CCGAGTCAATGAGATGAAAGCCGTGTGCTGTTGAAAGGCCAAGGCCCGGTGAAGGGCTCATCCACTTCTGTGAGAAAGCCAGGGCGGGGCGAG
```

Pick left primer, or use left primer below: Pick hybridization probe (internal oligo), or use oligo below: Pick right primer, or use right primer below (5' to 3' on opposite strand):

Pick Primers Reset Form

Sequence Id: A string to identify your output.

Targets: E.g. 50,2 requires primers to surround the 2 bases at positions 50 and 51. Or mark the [source sequence](#) with [and]: e.g. ...ATCT[CCCC]TCAT.. means that primers must flank the central CCCC.

Excluded Regions: E.g. 401,7 68,3 forbids selection of primers in the 7 bases starting at 401 and the 3 bases at 68. Or mark the [source sequence](#) with < and >: e.g. ...ATCT<CCCC>TCAT.. forbids primers in the central CCCC.

Product Size Ranges: 179-190

Number To Return: 5 Max 3' Stability: 9.0

Max Repeat Mispriming: 12.00 Pair Max Repeat Mispriming: 24.00

Max Template Mispriming: 12.00 Pair Max Template Mispriming: 24.00

Pick Primers Reset Form

General Primer Picking Conditions

Primer Size Min: 18 Opt: 20 Max: 27

Primer Tm Min: 65 Opt: 69 Max: 70 Max Tm Difference: 100.0 [Table of thermodynamic parameters:](#) Breslauer et al. 1986

Product Tm Min: Opt: Max:

Primer GC% Min: 20.0 Opt: Max: 80.0

Max Self Complementarity: 8.00 Max 3' Self Complementarity: 3.00

Max #N's: 0 Max Poly-X: 5

Inside Target Penalty: Outside Target Penalty: 0 [Note: you can set Inside Target Penalty to allow primers inside a target.](#)

First Base Index: 1 CG Clamp: 0

Concentration of monovalent cations: 50.0 Salt correction formula: Schildkraut and Lifson 1965

Concentration of divalent cations: 0.0 Concentration of dNTPs: 0.0

Annealing Oligo Concentration: 50.0 (Not the concentration of oligos in the reaction mix but of those annealing to template.)

Liberal Base Show Debugging Info Do not treat ambiguity codes in libraries as consensus Lowercase masking

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The screenshot shows the NCBI BLAST website. At the top, there is a navigation bar with the BLAST logo and links for Home, Recent Results, Saved Strategies, and Help. Below this, the main content area is divided into several sections:

- Basic Local Alignment Search Tool:** A section describing the tool's function: "BLAST finds regions of similarity between biological sequences. The program compares nucleotide or protein sequences to sequence databases and calculates the statistical significance." It includes a "Learn more" link.
- NEWS:** A box announcing a new version of IgBLAST (1.17) and a new field for "V frame shift" in the output.
- Web BLAST:** A section with three main options: "blastx" (translated nucleotide to protein), "tblastn" (protein to translated nucleotide), and "Protein BLAST" (protein to protein).
- BLAST Genomes:** A search box for selecting a genome, with options for Human, Mouse, Rat, and Microbes.
- Standalone and API BLAST:** Three options: "Download BLAST", "Use BLAST API", and "Use BLAST in the cloud".
- Specialized searches:** Four buttons for "SmartBLAST", "Primer-BLAST", "Global Align", and "CD-search", each with a brief description of its function.

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The image shows the NCBI Primer-BLAST web interface. The browser address bar shows the URL: https://www.ncbi.nlm.nih.gov/tools/primer-blast/index.cgi?LINK_LOC=BlastHome. The page title is "Primer-BLAST" and the subtitle is "A tool for finding specific primers". The main heading is "Finding primers specific to your PCR template (using Primer3 and BLAST)".

The interface is divided into several sections:

- PCR Template:** Includes a text input field for "Enter accession, gi, or FASTA sequence" and a "Range" section with "From" and "To" fields for "Forward primer" and "Reverse primer".
- Primer Parameters:** This section is highlighted with a red box. It contains:
 - "Use my own forward primer (5'->3' on plus strand)" with the input field containing "TCGGTTGCGTCGGCGTCT".
 - "Use my own reverse primer (5'->3' on minus strand)" with the input field containing "CCCGCCCTTGCCTTCTGC".
 - "PCR product size" with "Min" (70) and "Max" (1000) fields.
 - "# of primers to return" with a field containing "10".
 - "Primer melting temperatures (T_m)" with "Min" (57.0), "Opt" (60.0), "Max" (63.0), and "Max T_m difference" (3) fields.
- Exon/intron selection:** Includes a dropdown for "Exon junction span" (set to "No preference") and "Exon junction match" with "Min 5' match" (7), "Min 3' match" (4), and "Max 3' match" (8) fields. There is also an "Intron inclusion" checkbox and an "Intron length range" section with "Min" (1000) and "Max" (1000000) fields.

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The screenshot shows the NCBI Primer-BLAST web interface. The URL is https://www.ncbi.nlm.nih.gov/tools/primer-blast/index.cgi?LINK_LOC=BlastHome. The browser shows a OneDrive tab for 'Agrobiologie.cz - rozc...'. The interface includes several sections:

- Intron inclusion:** Minimal and maximal number of bases that must anneal to exons at the 5' or 3' side of the junction (7, 4, 8).
 Primer pair must be separated by at least one intron on the corresponding genomic DNA.
- Intron length range:** Min: 1000, Max: 1000000.
- Note:** Parameter values that differ from the default are highlighted in yellow.
- Primer Pair Specificity Checking Parameters:**
 - Enable search for primer pairs specific to the intended PCR template.
 - Search mode:** Automatic (highlighted in yellow).
 - Database:** Custom (highlighted in yellow).
Enter accession number, gi, or FASTA sequence (Clear):
NC_051806.1
NC_051807.1
NC_051808.1
NC_051809.1
NC_051810.1
NC_051811.1 (highlighted in red boxes)
 - Or, upload file: Procházet... Soubor nevybrán.
- Exclusion:**
 - Exclude predicted Refseq transcripts (accession with XM, XR prefix)
 - Exclude uncultured/environmental sample sequences
- Organism:** (highlighted in yellow)
Enter an organism name (or organism group name such as enterobacteriaceae, rodents), taxonomy id or select from the suggestion list as you type.
[Add more organisms](#)
- Entrez query (optional):** (highlighted in yellow)
- Primer specificity stringency:** Primer must have at least 2 total mismatches to unintended targets, including at least 2 mismatches within the last 5 bps at the 3' end. Ignore targets that have 6 or more mismatches to the primer.
- Max target amplicon size:** 1000 (highlighted in yellow)
- Allow splice variants:** Allow primer to amplify mRNA splice variants (requires refseq mRNA sequence as PCR template input)

At the bottom, there are two checked checkboxes: Show results in a new window and Use new graphic view. A final note states: **Note:** Parameter values that differ from the default are highlighted in yellow.

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The screenshot shows the NCBI Primer-BLAST web interface. The browser address bar displays the URL <https://www.ncbi.nlm.nih.gov/tools/primer-blast/primertool.cgi>. The page header includes the NIH logo and the text "U.S. National Library of Medicine" and "NCBI National Center for Biotechnology Information". The main heading is "Primer-BLAST" with the subtitle "A tool for finding specific primers". Below this, there is a sub-heading "Making primers specific to your PCR template. [more...](#)". A table shows the status of a submission:

Status	Submitted	Check
Current time	21 January 2021, 11:40:34	

At the bottom of the page, there is a footer with the text "BLAST is a registered trademark of the National Library of Medicine" and "NCBI National Center for Biotechnology Information, U.S. National Library of Medicine 8600 Rockville Pike, Bethesda MD, 20894 USA". There are also links for "Policies and Guidelines" and "Contact", and logos for "Support center", "Mailing list", "You Tube", "NATIONAL LIBRARY OF MEDICINE", "NIH", and "USA.gov".

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Primer-BLAST » JOB ID:gYtexPcg-ojdsmq3Z9dOhR3MX7cw30SqMQ

Primer-BLAST Results

Input PCR template: none
 Specificity of primers: Target templates were found in selected database: Custom
 Other reports: Search Summary

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GenBank number: [Browse GenBank]

or paste in your DNA sequence: (plain or FASTA format)

```
GCGCCGCTCGGTTTCGGTCGGCTGTCTGCGGCTCCTCCGCTGCCCTCGGACTCTGCAGCCCTCGGA
GCGCGCCCGGAGCGCGCCCGCCCTCGGAGTCATGGAGATGAAGCCGTGTGGTGTGGAAGGGCCAGGGC
CCGGTGGAGGGCACCATCCACTTCGTGCAGAGGCAAGGGCGGGGCGGAG
```

Submit

The sequence is: Linear Circular

Enzymes to use: NEB enzymes All commercially available specificities All specificities All + defined oligonucleotide sequences Only defined oligonucleotide sequences [define oligos]

Minimum ORF length to display: 100 a.a.

Name of sequence: WILD (optional)

Earlier projects:

Note: Your earlier projects will be deleted 2 days after they were last accessed. You need to have cookies enabled in your browser for this feature to work. Delete projects

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Linear Sequence: WILD

Display: - All commercial single cutter restriction enzymes
- Main non-overlapping, min. 100 aa ORFs

GC=73%, AT=27%

Legend:

Cleavage code	Enzyme name code
✂ blunt end cut	Available from NEB
⏏ 5' extension	Has other supplier
⏏ 3' extension	Not commercially available
⏏ cuts 1 strand	*: cleavage affected by CpG meth.
	#: cleavage affected by other meth.
	(enz.name): ambiguous site

Control Panel:

- Main options:** New DNA, Custom digest, View sequence, ORF summary, Save project, Print
- Availability:** NEB, All
- Display:** 2 cutters, 3 cutters
- Zoom:** Zoom in, More...
- List:** 0 cutters, **1 cutters** (highlighted), All sites, Save all sites, Flanking enzymes

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nc2.neb.com/NEBcutter2/listbycuts.php?name=0a69fde9-WILD&numcuts=1

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Single cutters [Help](#) [Comments](#)

WILD

Number of cuts = OK Sort order: Alphabetical Save as text file

#	Enzyme	Specificity	Sites & flanks	Cut positions (blunt - 5' ext. - 3' ext.)
1	AgsI	TT ₂ S [*] AA	list	132/131
2	Apal	G ₂ GGCC [*] C	list	*#145/141
3	BanI	G [*] GYRC ₂ C	list	154/158
4	BclI	CCATC [*] NNN [*] N ₂	list	166/167
5	BsgI	GTGCAG (N) ₁₄ NN [*]	list	190/188
6	BssHIII	G [*] CGCC ₂ C	list	*85/89
7	BstNI	CC [*] W ₂ GG	list	139/140
8	BtsCI	GGATG ₂ NN [*]	list	158/156
9	CviAII	C [*] AT ₂ G	list	105/107
10	DraIII	CAC ₂ NNN [*] GIG	list	*168/165
11	EcoO109I	RG [*] GNC ₂ CY	list	#141/144
12	EcoP15I	CAGCAG (N) ₂₅ NN ₂	list	28/30
13	FaiI	YA [*] TR	list	106
14	FatI	[*] CATG ₂	list	104/108
15	FokI	GGATG (N) ₉ NNNN ₂	list	145/149
16	Hinfi	G [*] ANT ₂ C	list	*101/104
17	Hpy188III	TC [*] NN ₂ GA	list	*98/100
18	MlyI	GAGTC (N) ₅ [*]	list	110
19	MmeI	TCCRAC (N) ₁₅ NN [*]	list	*48/46
20	MslI	CAYNN [*] NNRTG	list	109
21	MspAII	CMG [*] CKG	list	*44
22	NlaIII	₂ CATG [*]	list	108/104
23	NruI	TCG [*] CGA	list	*99
24	PleI	GAGTC [*] NNNN [*] N ₂	list	109/110
25	PspGI	[*] CCWGG ₂	list	#137/142
26	PspOMI	G [*] GGCC ₂ C	list	*#141/145
27	PstI	C ₂ TGCA [*] G	list	63/59
28	SfiI	C [*] TRYA ₂ G	list	59/63
29	TauI	G ₂ CSG [*] C	list	*35/32
30	TspDTI	ATGAA (N) ₉ NN [*]	list	127/125

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nc2.neb.com/NEBcutter2/listbycuts.php?name=0a69fde9-MUT&numcuts=1

Dokumenty - OneDrive Agrobiologie.cz - rozc...

Single cutters

MUT

Number of cuts = 1 OK Sort order: Alphabetical Save as text file

#	Enzyme	Specificity	Sites & flanks	Cut positions (blunt - 5' ext. - 3' ext.)
1	AgsI	TT_S^AA	list	132/131
2	Apal	G_GGCC^C	list	*#145/141
3	BaeGI	G^GCM^C	list	145/141
4	BccI	CCATCNNNN^N	list	166/167
5	BsgI	GTGCG(N)_14^NN^	list	190/188
6	BssHII	G^CGCG_C	list	*85/89
7	BstNI	CC^W^GG	list	139/140
8	BtsCI	GGATG_NN^	list	158/156
9	CviAII	C^AI_G	list	105/107
10	DraIII	CAC_NNN^GTG	list	*168/165
11	EcoO109I	RG^GNC_CY	list	#141/144
12	EcoP15I	CAGCAG(N)_25^NN^	list	28/30
13	FaiI	YA^TR	list	106
14	FatI	^CATG^	list	104/108
15	FokI	GGATG(N)_9^NNNN^	list	145/149
16	HinfI	G^ANT_C	list	*101/104
17	Hpy188III	TC^NN_GA	list	*98/100
18	MlyI	GAGTC(N)_5^	list	110
19	MmeI	TCCRAC(N)_18^NN^	list	*48/46
20	MslI	CAYNN^NNRTG	list	109
21	MspAII	CMG^CRG	list	*44
22	NlaIII	^CATG^	list	108/104
23	NruI	TCG^CGA	list	*99
24	PleI	GAGTCNNNN^N	list	109/110
25	PspGI	^CCWGG^	list	#137/142
26	PspOMI	G^GGCC_C	list	*#141/145
27	PstI	C^TGCA^G	list	63/59
28	SfiI	C^TRYA_G	list	59/63
29	TauI	G_CSG^C	list	*35/32
30	TspDTI	ATGAA(N)_9^NN^	list	127/125

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